

Medical Planning for Boldly Going Where No One Has Gone Before

National Aeronautics and
Space Administration



Benjamin Easter, MD, MBA, FACEP

Deputy Element Scientist
Exploration Medical Capability (ExMC)
NASA Human Research Program

Assistant Professor of Emergency Medicine
University of Colorado School of Medicine

Florida Atlantic University Emergency Medicine Grand Rounds
15 February 2022





Adapted Length-Based Resuscitation Tape - Adapted from Broselow Pediatric Emergency Tape

Gray	Pink	Red	Purple	Yellow	White	Blue	Orange
3-5 kg	Small Infant 6-7 kg	Infant 8-9 kg	Toddler 10-11 kg	Small Child 12-14 kg	Child 15-18 kg	Child 19-23 kg	Large Child 24-29 kg
50	50	50	60	60	60	70	80
1 straight	1 Straight	1 Straight	1 Straight	2 Straight	2 Straight	2 Straight or Curved	2 Straight or Curved
3.0 Cuffed	3.5 Uncuffed 3.0 Cuffed	3.5 Uncuffed 3.0 Cuffed	4.0 Uncuffed 3.5 Cuffed	4.5 Uncuffed 4.0 Cuffed	5.0 Uncuffed 4.5 Cuffed	5.5 Uncuffed 5.0 Cuffed	6.0 Cuffed
10-11	10-11	10-11	11-12	13-15	14-15	16-17	17-18
8	8	8	10	10	10	10	10
22-24	22-24	22-24	20-24	18-22	18-22	18-20	18-20
18/15	18/15	18/15	15	15	15	15	15
5-8	5-8	5-8	8-10	10	10	12-14	14-18
8	8	8	8-10	10-12	10-12	10-12	12
10-12	10-12	10-12	16-20	20-24	20-24	24-32	28-32

Guidelines for Endotracheal Tube Size and Distance:

Age	ET tube insertion distance (>2yr) = 3 x I
1 month	3
1 month to 8 months	3
8 months to 1 yr	4
1 yr & over	4
Tube Size Uncuffed	1
1 yr & older	1
Tube Size Cuffed	1



What's the problem space?

Now imagine if your ED was a Mars Medical System




- Pick everything you want for the next three years
- No more than 50 kg (notional)
- No resupply
- No real-time assistance (consults, UpToDate)
- No transfers



Agenda

- Designing Medical Systems for Human Spaceflight
- Martian Medical Analogue and Research Simulation
- How to get involved

Risk of Human Spaceflight



Medical events in and around spaceflight



Medical Condition	Events	Medical Condition	Events
Allergic reaction (mild to moderate)	11	Mouth ulcer	9
Ankle sprain/strain	11	Nasal congestion (space adaptation)	389
Back injury	31	Neck injury	9
Back pain (space adaptation)	382	Nose bleed (space adaptation)	6
Barotrauma (ear/sinus block)	31	Otitis externa	3
Choking/obstructed airway	3	Otitis media	3
Constipation (space adaptation)	113	Paresthesias	26
Diarrhea	33	Pharyngitis	11
Elbow sprain/strain	12	Respiratory infection	33
Eye abrasion (foreign body)	70	Shoulder sprain/strain	22
Eye chemical burn	6	Sinusitis	6
Eye infection	5	Skin abrasion	94
Finger dislocation	1	Skin infection	13
Fingernail delamination (EVA)	16	Skin laceration	1
Gastroenteritis	4	Skin rash	94
Headache (CO2 induced)	20	Smoke inhalation	3
Headache (late)	49	Space motion sickness (space adaptation)	325
Headache (space adaptation)	233	Urinary incontinence (space adaptation)	5
Hemorrhoids	2	Urinary retention (space adaptation) – female	5
Herpes Zoster reactivation (shingles)	1	Urinary retention (space adaptation) – male	4
Indigestion	6	Urinary tract infection – female	5
Influenza	1	Urinary tract infection – male	4
Insomnia (space adaptation)	299	Visual impairment/increased intracranial pressure (space adaptation)	15
Insomnia (late)	133	Wrist sprain/strain	5
Knee sprain/strain	7		

Five Hazards of Human Spaceflight




1. Distance from Earth
 2. Hostile/Closed Environments
 3. Isolation and Confinement
-
4. Radiation
 5. Altered Gravity Fields



Trail Map

Low Earth Orbit



400 km

Moon



~400,000 km

Mars



~225,000,000 km (average)

1000x

560x

560,000x



Trail Map (by relative distances from Denver)

Backyard



Rocky Mountains




South Pole




Progressive Earth Independence



- Real Time Communications
- Real Time Mission Control
- Strong Consumables Resupply
- Evacuation Capability (1.5 – 36 hrs)



- Near Real Time Communications
- Near Real Time Mission Control
- **Limited** Consumables Resupply
- Evacuation Capability (72 – 144 hrs)




- **No** Real Time Communications
- Delayed Mission Support
- **No** Consumables Resupply
- **No** Evacuation Capability





Increasing exposure to Hazards

Risk Modeling of Loss of Crew from Medical Events



1/90 is where the Space Shuttle TOTAL risk for LOC fell at the end of the program (Bagian, JAMA Neurology, 2019)

Exploration Medical System Challenges



Influencers

- Mission duration and objectives
- Distance from Earth
- Acceptable mission risk
- Time to definitive medical care
- Health and performance needs
- Level of training for on-board medical provider(s)
- Technology developments and medical advances

How do we solve this?


- Trade Space Analysis




The Backpack Problem



Mars Surface “Gear Ratios”







We can't fly everything that we want...

Mass

Data


Power

Time

Volume

Skills/training

The Backpack Solution



- Optimization problem
- Maximize value subject to constraints
- Question: what combination of blocks yields the most money but does not exceed the weight limit?

Examples of backpack solutions



FIGURE 4.2. Mercury medical kit containing items such as saline solution, bandages, stimulants, and decongestants (Photo courtesy of NASA)



What is a “trade space”?



We cannot take everything we want or need for a medical system due to resource constraints.

We cannot research every conceivable technology due to time and funding constraints!

Optimization/Prioritization is necessary!


- Rigorous exercise to identify the most *acceptable* solution amongst a set of alternative solutions
- Enables a traceable and repeatable process
- Evaluations will assess metrics for:
 - Risk: Probabilities of loss of crew and loss of mission
 - Resource characterizations: mass, power, volume, etc.



What questions are in the medical trade space?

- What if spacewalks are added to the mission parameters?
- What if the mission were extended?
- What if mass allocation lowers by 25%?
- Which imaging technology is most useful? Ultrasound or X-ray?
- What is the medical risk of reducing crew's water rations from 3L to 2L daily?
- If I could reduce the mass or volume of a piece of medical equipment by 50%, would it make the cut to fly in space?

Probabilistic Risk Assessment





SKIN

Burns secondary to Fire
Skin Abrasion
Skin Laceration

EYES

Acute Glaucoma
Eye Corneal Ulcer
Eye Infection
Retinal Detachment
Eye Abrasion
Eye Chemical Burn
Eye Penetration

EARS, NOSE, THROAT

Barotrauma (sinus block)
Nasal Congestion (SA)
Nosebleed (SA)
Acute Sinusitis
Hearing Loss
Otitis Externa
Otitis Media
Pharyngitis

DENTAL

Abscess
Caries
Exposed Pulp
Tooth Loss
Crown Loss
Filling Loss

CARDIOVASCULAR

Angina/Myocardial Infarction
Atrial Fibrillation / Atrial Flutter
Cardiogenic Shock secondary to Myocardial Infarction
Hypertension
Sudden Cardiac Arrest
Traumatic Hypovolemic Shock

GASTROINTESTINAL

Constipation (SA)
Abdominal Injury
Acute Cholecystitis
Acute Diverticulitis
Acute Pancreatitis
Appendicitis
Diarrhea
Gastroenteritis
Hemorrhoids
Indigestion
Small Bowel Obstruction

Pulmonary

Choking/Obstructed Airway
Respiratory Infection
Toxic Exposure: Ammonia
Smoke Inhalation
Chest Injury

*SA – Space Adaptation

NEUROLOGIC

Space Motion Sickness (SA)
Head Injury
Seizures
Headache
Stroke
Paresthesia
Headache (SA)
Neurogenic Shock
VIIP (SA)

MUSKULOSKELETAL

Back Pain (SA)
Abdominal Wall Hernia
Acute Arthritis
Back Injury
Ankle Sprain/Strain
Elbow Dislocation
Elbow Sprain/Strain
Finger Dislocation
Fingernail Delamination (EVA)
Hip Sprain/Strain
Hip/Proximal Femur Fracture
Knee Sprain/Strain
Lower Extremity Stress fracture
Lumbar Spine Fracture
Shoulder Dislocation
Shoulder Sprain/Strain
Acute Compartment Syndrome
Neck Injury
Wrist Sprain/Strain
Wrist Fracture

PSYCHIATRIC

Insomnia (Space Adaptation)
Late Insomnia
Anxiety
Behavioral Emergency
Depression

GENITOURINARY

Abnormal Uterine Bleeding
Acute Prostatitis
Nephrolithiasis
Urinary Incontinence (SA)
Urinary Retention (SA)
Vaginal Yeast Infection

INFECTION

Herpes Zoster (shingles)
Influenza
Mouth Ulcer
Sepsis
Skin Infection
Urinary Tract Infection

IMMUNE

Allergic Reaction
Anaphylaxis
Skin Rash
Medication Reaction


ENVIRONMENT

Acute Radiation Syndrome
Altitude Sickness
Decompression Sickness (EVA)
Headache (CO₂)

Examples of Output Visualizations for Clinicians




Top 10 Conditions Leading to LOCL




Run 1, med system


Top 10 Conditions Leading to LOCL



Run 2, resources removed

How much risk is from medical problems?





“...[The] assumption has been that risk of vehicle system malfunction far outweighs the risk of human system failure...NASA buys down the risk of failure of the human system through rigorous selection of individuals designed to minimize medical issues and optimize available capability in flight.”

NASA SP-2017-0633

How do you visualize the medical trade space?




System Characteristic	System 1- Run 1	System 2-Run 2
Mass (kg)	115.2	106.6
Prob of Loss of Crew	0.0057	0.0061
Prob of EVAC	0.0999	0.101
Crew Health Index	0.929	0.928
Requirements not met	0	4
Conditions not addressed	0	32






Minimizing Medical Risk is NOT the Goal



Notional


A Medical System Alone is Insufficient





Quick Hits on Other ExMC Focus Areas









Figure 1. Advanced RF and optical communications technologies combined with using the areostationary orbit offer 100-1000x greater data return from Mars and nearly continuous availability.

Table 1. High-Performance Mars-Earth Trunk Line Capability

	Frequency Band	Maturity	S/C Aperture	S/C Txmt Power	Ground Receiver	Data Rate (@ 2 AU)
Current State-of-the-Art (MRO)	X-band	Operational	3 m	100 W	34 m DSN BWG antenna	1 Mb/s
Next-Generation Trunk Line Options	Ka-band	TRL 6	3 m	200 W	34 m DSN BWG antenna	5 Mb/s
		TRL 3-4	5 m	1 kW	34 m DSN BWG antenna	70 Mb/s
	Optical (1550 nm)	TRL 6 (DSOC; to fly on 2023 Psyche Discovery Mission)	22 cm	4 W	5 m ground telescope	1 Mb/s
		TRL 3	50 cm	50 W	12 m ground telescope	100 Mb/s

We're not bringing an Intensive Care Unit but...



Notional

Medical Technology Demonstrations



Autonomous Medical Officer Support (AMOS)



- AMOS used to perform ultrasound of the bladder and kidneys in **autonomous fashion with no preflight training and no support from the ground**
- **1st spaceflight performance** of fully autonomous, untrained imaging
- Awarded the International Space Station Research and Development **2020 Compelling Results Award-Human Health**

Hemocue®



- Device for **point-of-care** analysis of white blood cell count and differential
- Validated performance in flight using control solutions and a fingerstick blood sample
- **1st real-time hematology performed in space**



Contents lists available at [ScienceDirect](https://www.sciencedirect.com)

Life Sciences in Space Research

journal homepage: www.elsevier.com/locate/issr




Spaceflight validation of technology for point-of-care monitoring of peripheral blood WBC and differential in astronauts during space missions





Exploration Lab Analysis

- Currently, no ability to perform real-time lab analysis in flight
- Required Assays
 - Metabolic Panel (components)
 - Blood Counts (components)
 - Urine
- Desired Assays
 - Cardiac, Blood Gases, Coags, ID
- Challenges
 - Reduce/eliminate waste
 - Fluidics in microgravity
 - Maintain industry standards



Imaging in the Exploration Environment



- Why imaging?
 - Of the 100 conditions that we scope medical system for:
 - 44 amenable to US
 - 27 amenable to XR
 - Possibly therapeutic (renal stone propulsion trial)
- Current Capability
- Challenges
 - Training
 - Obtaining high quality images (COMfORT, AR/VR, automated)
 - Image interpretation (store and forward, AI/ML)

Questions?



benjamin.easter@nasa.gov
@beasterMD

Trademark Statement



Trade names and trademarks are used in this report for identification only. Their usage does not constitute an official endorsement, either expressed or implied, by the National Aeronautics and Space Administration.